



Lojer

FOR EASY CARE

A professionals choice-
New Afia medical
examination tables
are made to perform.



Lojer Afia 4050

versatile and multifunctional examination tables for professionals

Afia 4040
Afia 4050
Afia 4140
Medical examination tables

Lojer Afia 4050, Afia 4040 and Afia 4140 medical examination tables



Afia 4050
Gynecological examination table



Afia 4140
Gynecological examination table



Afia 4040
General medical examination table

- 200 kg lifting capacity
- Safe and long lasting, sturdy frame and castors with central locking system. Power cable long enough to be connected to the wall socket from each side of the table.
- Top ergonomic usability; incl. electric height adjustment with round switch bar, which enables correct working position and creates safety while transferring the patient.
- Functionality and versatility are guaranteed by adequate set of standard equipment and a large set of additional accessories
- Modern design and fresh colours; lots of different colours for upholstery
- Can be cleaned easily; electric parts are splash proof, upholstery materials are resistant to cleaning- and disinfecting agents.

Models: Afia 4140, Afia 4040 and Afia 4050

Height adjustment	50-95 cm
Width	65 cm
Length	193 cm
Length of the back part	68 cm
Weight (without additional accessories)	110 kg
Castors	diam.125 mm
Safe working load (SWL)	200 kg

Standard equipment

- castors with central locking system
- electric adjustment of height with a foot pedal (round switch bar) and with hand control unit
- electric adjustment of the upper back region of 78 degrees with hand control unit
- potential equalization
- paper-roll holder
- heel holders and a gynecological dish in models 4050 and 4140
- gas spring assisted gynecological back section adjustment of 45 degrees and mechanical hip section adjustment of 15 degrees in model 4050,

Options (set in a factory)

- castors with separate locking
- special upholstery materials
- Trendelenburg (0-17 degrees)
- battery back-up

Additional accessories, which require an accessory bar and radial clamp

- additional accessory bar, length 40 cm, maximum 6 pcs/table (3 pcs/side)
- radial clamp for accessories bar
- adjustable arm rest
- IV-pole
- anesthesia frame
- various accessories, used at operating tables

Other accessories, that don't need accessories bars

- side rails with fasteners
- calf supports
- neck pillow



Guarantee 3 years

Standard colours for upholstery materials of Afia-medical examination tables



Mango 641-1060 Perle 641-0702 Rubin 641-0927 Silvergrey 641-1061 Sky 641-1053 Rauchblau 641-1059

